



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Inkubator cart BIO-RAD SCANGEL (Reconditioned)

Used, very good technical condition,

French production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 100 VA,

The card incubator is a microprocessor-controlled benchtop laboratory instrument designed to incubate (at a specific temperature for a specific time) Bio-Rad identification cards for antigen-antibody complex studies,

By incubating the identification cards, the optimum reaction temperature is achieved in the microtubes,

Designed for incubation at 37°C,

Standard programme: 15 minutes at 37°C ($\pm 2^{\circ}\text{C}$),

Programmable incubation time: 1 - 60 minutes,

Capacity: 24 ID cards placed on 2 racks,

Time and temperature are displayed on the LCD screen in English, French, German, Italian, Spanish or Portuguese,

Audible alarms for end of operation or failures,

Included:

2 12-place card racks,

Power cord,

User manual in English or French (PDF),

Chamber dimensions: 20.5 x 20.5 x 7 cm,

External dimensions: 34.5 x 32.5 x 20 cm,

Weight: 7.5 kg,

Cat 02,

Has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan

