



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



USG GE LOGIQ 3 (Reconditioned)

Used, very good technical condition,

U.S.A. production,

Power supply: 230 V,

Frequency: 50/60 Hz,

The Logiq 3 Pro is a high quality ultrasound scanner that is compact and economical to purchase and operate for use in hospital wards and private practice,

Fast algorithms for 3-D presentation,

Full compatibility with computers and computer networks (DICOM),

English language keyboard on the control panel,

Imaging modes:

B-mode,

M-mode,

Color Doppler,

Power Doppler,

Panoramic,

Virtual Convex,

Harmonic,

Included:

3C LINEAR head,

GE E8C endovaginal transducer,

15" colour monitor,

SONY CDR-W Recorder,

Thermal printer - Videoprinter Mitsubishi P93,

Black and white,

Paper roll width 110 mm,

Two active head ports,

Transducer and gel holders,

Mobile, on castors with brakes,

Dimensions: 47 x 138 x 94 cm,

Weight of set: 162 kg,

Has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): instalments, leasing, loan

