



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Oxygen Therapy Kit 5L (Reconditioned)**

Cylinder capacity: 5 litres,  
Working pressure PW: 200 bar,  
Test pressure PH: 300 bar,  
Equipped with 3/4" valve,  
Certified and validated,  
Including a new oxygen regulator with flow control and humidifier (type of regulator available depending on stock) and a choice of:  
Oxygen mask for adults with 2.1 m tube,  
Paediatric oxygen mask with 2.1 m tubing,  
Nasal cannula for adults with 2.1 m tube,  
Cylinder dimensions:  
Height: 54.5 cm,  
Diameter: 13.5 cm,  
Pictures for reference,  
Comes with a current Technical Passport issued,  
Warranty:  
6 months for domestic market - Poland,  
1 month for international market,  
Financing options (Poland only): Installments, Leasing, Loan

