



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Linear head GE 7L (Reconditioned)

Used, good technical condition,

U.S.A. production,

Probe type: linear,

Frequency range: 2 - 7 MHz,

The GE 7L probe is an ultrasound transducer that is compatible with many GE ultrasound systems:

Vivid 7 Dimension,

Vivid 7 Pro,

Logiq 5 Expert,

Logiq 7,

Logiq E9,

Logiq S6,

Application: paediatric, small parts and vascular,

Dimensions:

Handle: 11 x 3 x 6.5 cm,

Cable length: 202 cm,

Weight: approx. 1.4 kg,

Cat 02,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan

