



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE T739 intraoperative head (Reconditioned)

Used, good technical condition,

U.S.A. production,

Type of probe: linear array,

Frequency range: 4 - 10 MHz,

The GE T739 probe is an ultrasound transducer that is compatible with many GE ultrasound systems:

Logiq 5 Expert,

Logiq 7,

Logiq P5,

Logiq P6,

Logiq E9,

Applications: vascular, small parts and intraoperative ultrasound images of both adults and children,

Dimensions:

Handle: 4.2 x 4.7 x 1.7 cm,

Cable length: 350 cm,

Weight: approx. 1.3 kg,

Comes with a current Technical Passport issued,

Guarantee:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan

