



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE 7S sectoral head (Reconditioned)

Used, very good technical condition,
U.S.A. production (2010),
Probe type: sector,
Frequency range: 3.7 - 6.3 MHz,
Field of view 90°,
Compatible with many GE ultrasound systems: are
Logiq 3,
Logiq 5 Expert,
Logiq 7,
Logiq P5,
Logiq P6 and others.
Specialising in cardiac and abdominal applications,
This transducer can produce cardiac images of both adults and children,
Dimensions:
Handle: 3.5 x 10.2 x 2.9 cm,
Cable length: 208 cm,
Weight: approx. 1.1 kg,
Has a current Technical Passport issued,
Guarantee:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan

