



Docking station for 4 pumps B.BRAUN 8713140 (Reconditioned)

Used, very good technical condition,

German production,

Power supply: 230 V,

Emergency power supply: built-in NiMH 4.8V / 1.9Ah battery - can be replaced at additional charge,

Frequency: 50/60 Hz,

Power: 80 VA,

It has a comfortable handle for movement,

The station is equipped with an acoustic alarm system,

The station is equipped with an acoustic alarm system, the alarm volume is adjustable,

The station is equipped with an acoustic alarm system, The station has a volume control for alarms, It has a mounting location for the bracket with which it can be mounted on a tripod (no bracket included),

No pumps included - can be retrofitted with infusion pumps (volumetric) B. BRAUN INFUSOMAT SPACE for an additional charge,

Dimensions: 15 x 27 x 34 cm.

Weight: 4 kg,

Cat 03,

Has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing: Installments, Leasing, Loan.

