



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Camera Processor OLYMPUS OTV-SE (Reconditioned)

- Used, very good technical condition,
- Japanese production,
- Power supply 230 V,
- Frequency 50-60 Hz,
- Dimensions: 22 x 29,5 x 6,5 cm,
- It has exoised the current Technical Passport,
- Warranty:
 - 6 months for national field - Poland,
 - 1 month for international field,
- The possibility of financing (Polish only): Installments, lease, loan.

