



Laparoscopy pump AESCULAP PROMIS LINE FLUID 2 (Reconditioned)

Used, very good technical condition,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power consumption: 70 W,

PROMIS LINE FLUID 2 is an advanced system ensuring maximum patient safety. Thanks to its simple and intuitive operation, it improves comfort and safety at work, thus ensuring full functionality,

Maximum flushing and suction power: 500 mmHg,

Flow range: 0 - 1000 ml / min,

Instruction manual in English (PDF),

Dimensions: 32 x 30 x 13 cm,

Weight of set: 9 kg,

Cat 02,

Has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan

