



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS endoscope kit (CLV-S40, OTV-S7) (Reconditioned)

Used, very good technical condition,

Japanese production,

The set includes:

OLYMPUS VISERA OTV-S7 Image processor:

Used, perfect technical condition,

Japanese production,

Power supply 230 V,

Frequency 50/60 Hz,

Power: 70 VA,

Keypad on control panel in English,

Foiled control panel,

TV type: PAL,

Endoscope functions:

WHITE BALANCE function - Automatic white balance,

+/- 8 steps light intensity adjustment,

Freeze function,

Outputs:

1 x BNC output,

2 x Y/C output,

RGB,

H.COPY,

VTR REMOTE,

SYSTEM,

KEYBOARD,

LIGHT CONTROL,

Adjustable focus,

Dimensions: 29,5 x 38 x 9,5 cm,

It has a current Technical Passport,

OLYMPUS Light Source Model VISERA CLV-S40:

Used, very good technical condition,

Japanese production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 3 A,

Keypad on control panel in English,

Foiled control panel,

Suitable for flexible and rigid endoscopes,

Xenon light,

Bulbs:

Xenon MD-631 of 300 W - main light,

Halogen 150 W - emergency lighting,

Service life of the lamp: approx. 500 hours of continuous use,

Built-in air pump with electronic power control,

Automatic brightness control,

16-stage light intensity control,

Outputs:

VIDEO IN,

VIDEO OUT,

H.COPY,

SYSTEM,

LIGHT CONTROL,

Dimensions: 37 x 41 x 14 cm,

It has a current Technical Passport issued,

OLYMPUS Keyboard Model N860-8745-T402, MAJ-1124 for endoscopic camera OTV-S7V,

Light Control Cable (MH-988),

User manual in English (PDF),

Set weight: 19.5 kg,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing option (Poland only): instalments, leasing, loan

