



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Headrest FAMED ŻYWIEC WS-58.5 (Reconditioned)**

Used, very good technical condition,  
Made in Poland,  
Mattress made of polyurethane in black colour,  
Metal construction made of stainless steel,  
Integrated fixing block,  
Mattress surface dimensions: 52 x 29 x 5 cm,  
The mattress has a small hole at the bottom of the upholstery, the damage does not affect the use in any way,  
Dimensions: 52 x 41.5 x 10 cm  
Weight: 7 kg,  
Has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan

