



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VOLKER 3080 S Rehabilitation bed (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Germany,
- Power supply: 230 V,
- Power: 150 VA,
- Völker 3080 S is a hospital bed that meets a wide range of requirements for common use in clinics and hospitals,
- With the successful combination of form and function, the bed can provide greater comfort to the patient while allowing caregivers to work safely and efficiently,
- Allows patients to go easier through the various steps of recovery, and to the overall better well-being of the patient. With the support of caregivers - from sitting in the Comfort Chair position through getting up and out of bed safely,
- Certified to IPX4 protection category (waterproof measurement standard),
- Patented Völker side rail system:
 - Four independent side rails that move according to the profile of the lying surface, which are divided on both sides and can be adjusted to one or two levels,
 - The main advantage of this solution is that users are always protected from falling from the bed, even in a sitting position. As a result, patients feel adequately protected from accidentally rolling out of bed. The side rails are very versatile and durable, and can also help with rehabilitation,
- Völker MiS comfort microstimulation system:
 - The system detects even small movements made by the patient. The lying surface then reflects them back through the mattress to the patient's body through tactile stimulation. This responsive feedback supports spontaneous movements and microcirculation in the patient's tissues to improve blood circulation. This can lead to greater comfort, better sleep and increased mobility. Patients can sleep well and wake up rested,
- Features:
 - 3-segment,
 - Remote control and pull-out control panel for staff,
 - It is possible to lock the functions on the remote control with a locking hook attached to the remote control cable,
 - Electrically controlled: up/down, backrest, foot segment, anti-Trendelenburg position,
 - Additional manual control of backrest and foot segment,
 - Electric bed height adjustment in the range: 42 - 80 cm,
 - The backrest can be raised to a maximum angle of 70°,
 - The upper part of the lying surface leg can be raised to an angle of 60°,
 - Trendelenburg and anti-Trendelenburg position in the range up to: +/- 30°,
 - Has 4 built-in slots for IV extension arms and patient handles,
- Mobile on 4 retractable wheels,
- User manual in English (PDF),
- Bed dimensions: 188 x 85 x 1.5 cm,
- Overall dimensions: 210 x 98 x 78 / 116 cm,
- Maximum load: 170 kg,
- Weight: 94 kg,
- It has a current Technical Passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,
- Valid technical passport issued,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

