



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Docking station for 6 pumps FRESENIUS KABI ORCHESTRA BASE A (Reconditioned)

Used, very good technical condition,
French production,
Power supply: 230 V
Emergency power supply: through built-in 6 V batteries - replaced with new ones,
Frequency: 50/60 Hz,
Power: 125 VA,
The Base A docking station is a data hub and power supply for up to 6 DPS modules and/or DPS Visio modules and/or MVP modules,
The station has a crack on the housing, the crack does not affect the use of the device in any way,
Allows faster therapeutic control,
Communication interface compatible with PDMS main software,
Built-in RS232 interface,
One single power cable for all connected modules,
Minimum battery life of 2 hours,
Equipped with tripod with comfortable carrying handle,
No pumps included - can be retrofitted to FRESENIUS ORCHESTRA infusion pumps (syringe or drip) at additional cost
Operating instructions in German (PDF),
Docking station dimensions:
Height: 8 cm,
Width: 34.5 cm,
Depth: 25.5 mm,
Height including tripod: 80 cm,
Weight of station with tripod: 4 kg,
Weight of docking station with pumps: 17.5 kg,
Cat 03,
It has an up-to-date Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan

