



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



OLYMPUS CV-145 video processor with CLE-145 light source (set) (Reconditioned)

Used, very good technical condition,

Japanese production,

The set includes:

OLYMPUS Image Processor Model CV-145 EVIS EXERA:

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 0.5 A,

Automatic white balance control,

It has a "freeze" function,

Colour saturation adjustment: red and blue,

Foiled control panel,

Dimensions: 42 x 45 x 7.5 cm,

Weight: 7.8 kg,

Light source OLYMPUS Model CLE-145 EVIS EXERA:

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 1.5 VA,

Suitable for flexible endoscopes,

Halogen bulb: 15V 150W,

Emergency lighting: halogen 15V 150W,

Foiled, touch control panel,

Automatic/manual regulation of image brightness,

Built-in air pump with electronic power control,

Dimensions: 37 x 45 x 14.5 cm,

Weight: 12.8 kg,

Keyboard OLYMPUS MAJ - 1124,

Light Control cable MH - 966,

Bottle for rinsing the endoscope,

The set can be retrofitted with Pigtail MAJ - 843 for an additional charge,

Weight of set: 21.9 kg,

Cat 04,

It has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan

