



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Veterinary table VET-E01 (New)

New,
Made in Poland,
Veterinary treatment table model VET-E01, has a stable steel construction on adjustable anti-slip feet, which allows easy access to the table top,
The entire structure is coated with powder paint resistant to mechanical damage and scratches,
As standard, the table has a fully bearing-mounted electric height adjustment with a manual remote control and a table top covered with non-slip material. As a standard, the table is equipped with holders enabling the assembly of additional options,
Additional options are available (at an extra charge);
IV hanger,
Euro rails,
Tool table,
Foot remote control,
Longitudinal outrigger,
Cross boom,
Castors,
Tilting bowl,
Technical data:
Overall length: 135 cm,
Height: 56 - 106 cm,
Total width: 65 cm,
Table top type: Covered with anti-slip material,
Electric adjustment of the tilt angle of the top: No,
Maximum pressure: 150 kg,
Weight of device: 53 kg,
Electrical voltage: 230 V / 50-60 Hz,
Maximum power: 160 W,
Protection class: II,
Enclosure: IP54,
Warranty on structure 24 months,
Warranty for hydraulic accessories 12 months

