



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Veterinary table VET-E05 (New)

New,
Made in Poland,
Veterinary treatment table model VET-E05, has a stable steel construction on adjustable anti-slip feet, which allows easy access to the table top,
The entire structure is coated with powder paint resistant to mechanical damage and scratches,
As standard, the table is equipped with a fully bearing-mounted electric height adjustment with a manual remote control and a table top made of acid-resistant steel,
The VET-E05 model also enables adjustment of the table top inclination angle with an electric actuator. As a standard, the table is equipped with holders enabling the assembly of additional options.
Additional options are available (at an extra charge);
IV hanger,
Euro rails,
Tool table,
Foot remote control,
Longitudinal outrigger,
Cross boom,
Castors,
Mattress made of waterproof material,
Tiltable bowl,
Technical data:
Overall length: 135 cm,
Height: 56 - 106 cm,
Total width: 65 cm,
Type of top: made of acid-resistant steel,
Electric adjustment of the top's angle: Yes,
Range of top adjustment: 0°-30°,
Maximum pressure: 150 kg,
Weight of the machine: 53 kg,
Electrical voltage: 230 V / 50-60 Hz,
Maximum power: 160 W,
Protection class: II,
Enclosure: IP54,
Warranty on structure 24 months,
Warranty for hydraulic accessories 12 months

