



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Veterinary table VET-E06 (New)**

New,  
Made in Poland,  
Veterinary treatment table model VET-E06, has a stable steel construction on adjustable anti-slip feet, which enables easy access to the tabletop,  
The entire structure is coated with powder paint resistant to mechanical damage and scratches,  
As standard, the table is equipped with a fully bearing-mounted electric height adjustment with a manual remote control and a deepened table top made of acid-resistant steel,  
Additionally, in VET-E06 model, the table top inclination angle may be adjusted by means of an electric actuator. As a standard, the table is equipped with holders enabling the assembly of additional options.  
Additional options are available (at an extra charge);  
IV hanger,  
Euro rails,  
Tool table,  
Foot remote control,  
Longitudinal outrigger,  
Cross boom,  
Castors,  
Mattress made of waterproof material,  
Tiltable bowl,  
Technical data:  
Overall length: 135 cm,  
Height: 56 - 106 cm,  
Total width: 65 cm,  
Table top type: Deepened top made of acid-resistant steel,  
Electrical adjustment of the top angle: Yes,  
Range of top adjustment: 0°-30°,  
Maximum pressure: 150 kg,  
Weight of the machine: 53 kg,  
Electrical voltage: 230 V / 50-60 Hz,  
Maximum power: 160 W,  
Protection class: II,  
Enclosure: IP54,  
Warranty on structure 24 months,  
Warranty for hydraulic accessories 12 months

