



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Veterinary table VET-H01 (New)

New,

Made in Poland,

Veterinary treatment table model VET-H01, has a stable steel construction on adjustable anti-slip feet, which enables easy access to the tabletop,

The entire structure is coated with powder paint resistant to mechanical damage and scratches,

As a standard, the table has a fully bearing-mounted hydraulic height adjustment with a foot lever and a table top covered with an anti-slip material.

As standard, the table is equipped with handles enabling the assembly of additional options.

Additional options are available (at an extra charge);

IV hanger,

Euro rails,

Tool table,

Longitudinal outrigger,

Cross boom,

Castors,

Tilting bowl,

Technical data:

Overall length: 135 cm,

Height: 56 - 106 cm,

Total width: 65 cm,

Type of top: covered with non-slip material,

Manual adjustment of the tilt angle of the top: No,

Maximum pressure: 200 kg,

Weight of device: 53 kg,

Warranty on structure 24 months,

Warranty for hydraulic accessories 12 months

