



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Veterinary table VET-H03 (New)**

New,

Made in Poland,

Veterinary treatment table model VET-H03, has a stable steel construction on adjustable anti-slip feet, which enables easy access to the tabletop,

The entire structure is coated with powder paint resistant to mechanical damage and scratches,

As a standard, the table has a fully bearing-mounted hydraulic height adjustment with a foot lever and a deepened top made of acid-resistant steel.

As a standard, the table is equipped with holders enabling the assembly of additional options.

Additional options are available (at an extra charge)

IV hanger,

Euro rails,

Tool table,

Longitudinal outrigger,

Cross boom,

Castors,

Mattress made of waterproof material,

Tilting bowl,

Technical data:

Overall length: 135 cm,

Height: 56 - 106 cm,

Total width: 65 cm,

Type of top: Deepened top made of acid-resistant steel,

Manual adjustment of the tilt angle of the top: No,

Maximum pressure: 200 kg,

Weight of device: 53 kg,

Warranty on structure 24 months,

Warranty for hydraulic accessories 12 months

