



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Veterinary table VET-H04 (New)

New,
Made in Poland,
Veterinary treatment table model VET-H04, has a stable steel construction on adjustable anti-slip feet, which enables easy access to the table top,
The entire structure is coated with powder paint resistant to mechanical damage and scratches,
As a standard, the table is equipped with a fully bearing hydraulic height adjustment with a foot lever and a table top covered with an anti-slip material,
The VET-H04 model also features gas spring adjustment of the table top angle.
As standard, the table is equipped with holders enabling the assembly of additional options,
Additional options are available (at an extra charge);
IV hanger,
Euro rails,
Tool table,
Longitudinal outrigger,
Cross boom,
Castors,
Tilting bowl,
Technical data:
Overall length: 135 cm,
Height: 56 - 106 cm,
Total width: 65 cm,
Type of top: covered with non-slip material,
Manual adjustment of the tilt angle of the top: Yes,
Top adjustment range: 0°-30°,
Maximum pressure: 200 kg,
Weight of the machine: 53 kg,
Warranty on structure 24 months,
Warranty for hydraulic accessories 12 months

