



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Veterinary table VET-H05 (New)

New,
Made in Poland,
Veterinary treatment table model VET-H05, has a stable steel construction on adjustable anti-slip feet, which enables easy access to the tabletop,
The entire structure is coated with powder paint resistant to mechanical damage and scratches,
As standard, the table is equipped with a fully bearing hydraulic height adjustment with a foot lever and a table top made of acid-resistant steel,
In addition, the VET-H05 model enables adjustment of the table top inclination angle with a gas spring actuator.
As a standard, the table is equipped with holders enabling the assembly of additional options,
Additional options are available (at an extra charge);
IV hanger,
Euro rails,
Tool table,
Longitudinal outrigger,
Cross boom,
Castors,
Mattress made of waterproof material,
Tilting bowl,
Technical data:
Overall length: 135 cm,
Height: 54 - 104 cm,
Total width: 65 cm,
Type of top: made of acid-resistant steel,
Manual adjustment of the top angle: Yes,
Range of top adjustment: 0°-30°,
Maximum pressure: 200 kg,
Weight of the machine: 53 kg,
Warranty on structure 24 months,
Warranty for hydraulic accessories 12 months

