



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **VET-RE1 veterinary X-ray table (New)**

New,  
Made in Poland,  
Veterinary treatment table VET-RE1, has a stable steel construction on adjustable anti-slip feet, which enables free access to the tabletop,  
The entire structure is coated with powder paint resistant to mechanical damage and scratches,  
The entire structure is covered with powder coating resistant to mechanical damage and scratches. As a standard, the table is fitted with ball-bearing electric height adjustment with a manual remote control and the top made of HPL material transparent to X-rays,  
As a standard, the table is fitted with handles enabling the assembly of additional options and trays for X-ray cassettes, which can be set along the entire length of the tabletop,  
As additional options are available (for an extra charge);  
IV hanger,  
Euro rails,  
Instrument table,  
Foot remote control,  
Longitudinal outrigger,  
Cross boom,  
Castors,  
Mattress made of waterproof material,  
Technical data:  
Overall length: 135 cm,  
Height: 58- 108 cm,  
Total width: 65 cm,  
Type of top: made of hpl plastic,  
Electric adjustment of the tilt angle of the top: No,  
Maximum pressure: 150 kg,  
Weight of device: 53 kg,  
Electrical voltage: 230 V / 50-60 Hz,  
Maximum power: 160 W,  
Protection class: II,  
Enclosure: IP54,  
Warranty on structure 24 months,  
Warranty for hydraulic accessories 12 months

