



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VET-RE2 veterinary X-ray table (New)

New,
Made in Poland,
Veterinary treatment table VET-RE1, has a stable steel construction on adjustable anti-slip feet, which enables free access to the tabletop, The entire structure is coated with powder paint resistant to mechanical damage and scratches,
As a standard, the table has a fully bearing-mounted electric height adjustment with a manual remote control and a tabletop made of HPL material transparent for X-rays,
Additionally, in VET-RE2 model, it is possible to adjust the table top inclination angle with an electric actuator.
As a standard the table is equipped with handles enabling assembly of additional options and trays for X-ray cassettes, which can be set along the entire length of the tabletop,
Additional options are available (at an extra charge);
IV hanger,
Euro rails,
Instrument table,
Foot remote control,
Longitudinal outrigger,
Cross boom,
Castors,
Mattress made of waterproof material,
Technical data:
Overall length: 135 cm,
Height: 58- 108 cm,
Total width: 65 cm,
Type of top: made of hpl plastic,
Electric adjustment of the top angle: Yes,
Top adjustment range: 0° - 30°,
Maximum pressure: 150 kg,
Weight of the machine: 53 kg,
Electrical voltage: 230 V / 50-60 Hz,
Maximum power: 160 W,
Protection class: II,
Enclosure: IP54,
Warranty on structure 24 months,
Warranty for hydraulic accessories 12 months

