



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VET-RH1 veterinary X-ray table (New)

New,

Made in Poland,

Veterinary treatment table VET-RH1, has a stable steel construction on adjustable anti-slip feet, which allows easy access to the table top, The entire structure is coated with powder paint resistant to mechanical damage and scratches,

The entire structure is covered with powder coating resistant to mechanical damage and scratches. As a standard, the table is fitted with a hydraulic height adjustment with a foot lever and the top made of HPL material transparent to X-rays,

As a standard, the table is equipped with holders enabling the assembly of additional options and trays for X-ray cassettes, which can be set along the entire length of the tabletop,

As additional options are available (for an extra charge);

Hanger for an IV,

Euro rails,

Instrument table,

Longitudinal outrigger,

Cross boom,

Castors,

Mattress made of waterproof material,

Technical data:

Overall length: 135 cm,

Height: 58- 108 cm,

Total width: 65 cm,

Type of top: made of hpl plastic,

Manual adjustment of the tilt angle of the top: No,

Maximum pressure: 150 kg,

Weight of device: 53 kg,

Warranty on structure 24 months,

Warranty on hydraulic accessories 12 months

