



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VET-RH2 veterinary X-ray table (New)

New,
Made in Poland,
Veterinary treatment table VET-RH2, has a stable steel construction on adjustable anti-slip feet, which allows easy access to the table top,
The entire structure is coated with powder paint resistant to mechanical damage and scratches,
As a standard, the table has a fully bearing-mounted hydraulic height adjustment with a foot lever. The table has a fully bearing-mounted hydraulic height adjustment with a foot lever and the top is made of HPL material transparent to X-rays,
Additionally, in VET-RH2 model the tabletop angle can be adjusted with a gas spring,
As a standard, the table is equipped with handles enabling assembly of additional options and trays for X-ray cassettes, which can be set along the entire top length,
As additional options are available (for an extra charge);
IV hanger,
Euro rails,
Instrument table,
Longitudinal outrigger,
Cross boom,
Castors,
Mattress made of waterproof material,
Technical data:
Overall length: 135 cm,
Height: 58- 108 cm,
Total width: 65 cm,
Type of top: made of hpl plastic,
Manual adjustment of the top angle: Yes,
Top adjustment range: 0° - 20°,
Maximum pressure: 150 kg,
Weight of the device: 53 kg,
Warranty on structure 24 months,
Warranty on hydraulic accessories 12 months

