



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VET-RS1 veterinary X-ray table (New)

New,
Made in Poland,
Veterinary treatment table VET-RS1, has a stable steel construction on adjustable anti-slip feet, which allows easy access to the table top,
The entire structure is coated with powder paint resistant to mechanical damage and scratches,
The table has a fixed height and a top made of HPL material transparent for X-rays.
Additional options are available (at an extra charge)
Hanger for an IV,
Euro rails,
Tool table,
Longitudinal outrigger,
Cross boom,
Castors,
Mattress made of waterproof material,
Technical data:
Overall length: 135 cm,
Height: 92 cm,
Total width: 65 cm,
Type of top: made of hpl plastic,
Manual adjustment of the tilt angle of the top: No,
Maximum pressure: 150 kg,
Weight of device: 53 kg,
Warranty on structure 24 months,
Warranty on hydraulic accessories 12 months

