



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment couch VET-L03 (New)

New,
Made in Poland,
Veterinary treatment couch model VET-L03 has a stable, steel construction of fixed height with possibility of levelling by means of unscrewable foot,
The entire structure is coated with powder paint resistant to mechanical damage and scratches,
The couch has a deepened table top made of acid-resistant steel,
Additional options are available (at an extra charge);
IV hanger,
Euro rails,
Tool table,
Longitudinal outrigger,
Cross boom,
Castors,
Mattress made of waterproof material,
Technical data:
Overall length: 135 cm,
Height: 85 cm,
Total width: 65 cm,
Type of top: Deepened top made of acid-resistant steel,
Manual adjustment of the tilt angle of the top: No,
Maximum pressure: 150 kg,
Weight of device: 53 kg,
Warranty on structure 24 months,
Warranty on hydraulic accessories 12 months

