



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Veterinary stretcher (New)

New,
Made in Poland,
Veterinary stretcher - Stretcher has a steel construction covered with powder paint resistant to mechanical damage and scratches, filling from scratch-resistant and easily washable material,
Technical data:
Overall length: 122 cm,
Total width: 62 cm,
Manual adjustment of the tilt angle of the top: No,
Maximum pressure: 150 kg,
Unit weight: 53 kg,
Warranty: 24 months

