



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PENTAX KW2422S biopsy forceps (Reconditioned)

Used, very good technical condition,
Japanese production,
Working channel diameter: 2,8 mm,
Working length: 2200 mm,
Oval spoon with window,
Without needle,
Autoclavable,
Visible signs of use on the handle,
Cat 02

