



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDTRONIC (PHYSIO-CONTROL) LIFEPAK 15 defibrillator bed connector (Ex-Display)

New, demo,
Perfect technical condition,
U.S.A. production,
The bed connector attaches the LIFEPAK 15 defibrillator/monitor to a stretcher or bed while transporting patients inside the hospital or on the way from the ambulance to the hospital,
The connector is recommended for hospital use only,
Dimensions: 21.5 x 15 x 6.5 cm,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Instalments, Leasing, Loan

