



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FRESENIUS ORCHESTRA MODULE MVP PT DIN Infusion Pump (volumetric) (Reconditioned)

Used, very good technical condition,
French production,
Power supply: 230 V,
Emergency power supply: through a built-in 6 V battery - can be replaced at extra cost,
Frequency: 50/60 Hz,
Power: 10 W,
Orchestra MVP PT module is a volumetric infusion pump that can be used with all Orchestra configurations,
Foiled control panel,
Selectable menu languages: PL, EN, FR, DE, IT, SP, NL, PT, SV, NOV,
Keypad on control panel in English,
Flow rate range from 1 to 1000 ml / h in normal mode, from 0.1 to 100 ml / h in micro-infusion mode
Flow rate accuracy $\pm 5\%$ with recommended sets,
Volume range from 1 to 9999 ml in normal mode, 0.1 to 999.9 ml in micro-infusion mode,
Infusion time adjustable from 1 minute to 96 hours,
Infusion modes: micro or macro infusion, increase/decrease, sequential, induction, bolus, primary/secondary,
DPS (Dynamic Pressure System): detection of fluctuations in line pressure, allowing both faster detection of occlusion and line disconnection,
Anti-Bolus System: automatic bolus reduction when occlusion is released,
Pressure limit adjustable from 100 to 900 mmHg,
Battery life: minimum 2 hours at 125 ml / h,
Audible alarms,
Equipment:
AC adapter,
Droplet sensor,
Instruction manual in German (PDF),
Dimensions: 35 x 12.5 x 13.5 cm,
Weight: 2.5 kg,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan

