



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Certified Test / COVID Tests - 19 IgM/IgG LANSIONBIO (25 pcs / pack) (New)

New,
Manufactured by P.R.C.,
CE approved,
For professional use only,
This test kit is designed for the in vitro diagnosis of COVID-19 IgM/IgG antibodies in human serum, plasma and whole blood using the immunofluorescence method,
When a test strip is inserted into the analyzer, the instrument automatically scans the ribbons and detects the fluorescence intensity of the composite emission from the test and control areas. Analysis of the fluorescence value is used to calculate the content of the detected substances,

Specifications:

Response time: 15 minutes

Sample: serum, plasma and whole blood,

Methodology: Fluorescence immunoassay,

Storage : 4-30 ° C,

Included:

IgM / IgG test strip in sealed bag with desiccant - 25 pcs.

Sample diluent - 1 pc.

User Manual in English - 1 item.

Card with QR code for calibration - 1 item.

Certificate,

The tests are adapted for LANSIONBIO analyzers:

LS - 1000,

LS - 2000,

LS - 1100,

LS - 2100,

LS - 4000,

LS - 7000,

Insert the test strip into the analyzer immediately after taking the sample. The result will be displayed on the screen and automatically printed.

