



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Connecting cable for electrosurgery ERBE 20196-65 (New)

New,
German production (2019),
ERBE connecting cable transfers electricity from equipment to instruments without changing the energy,
Designed exclusively for Erbe equipment,
Specifications:
Bipolar cable,
8/4 mm,
For Aesculap tweezers,
Bent connector,
Cable length: 4.5 m,
Maximum peak voltage: 1 kVp
Basic Features:
Compatibility : VIO, ICC, ACC, T series,
Connection : standard (BI 8/4),
Cleaning and preparation for use:
Machine cleaning: max 95 °C,
Autoclave sterilization: max 138 °C,
Multilingual instruction manual,
Warranty: 12 months

