



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **CLINNOVA Clinical 2 examination and treatment table (NV2560-PRM-SEERSMEDICAL) (New)**

New,  
English production,  
Very strong and rigid table available with a wide range of accessories,  
Suitable for very heavy patients during examinations and minor medical procedures,  
Height adjustment of the table with an electric actuator, manual adjustment of the backrest with a gas spring,  
The table is equipped with a premium base and paper holder as standard,  
Table advantages:

Retractable wheel design and adjustable leg for stability,  
Frame construction designed for easy use of mobile lift,  
Antimicrobial powder-coated base frame,  
Easy Clean 70 cm wide upholstery as standard,  
Choice of 15 standard upholstery colours  
The table rises from minimum to maximum height in only 15 seconds,  
Handheld remote control fitted as standard with this model,

### Technical data:

Safe working load and lifting capacity of 265 kg,  
Adjustable backrest from -20° to 80° with anti-trap design,  
Standard width of the bed is 70cm made of foam,  
Table height adjustment range from 46 cm to 95 cm,

### Additional accessories (at extra charge):

Face hole with plug in the bolster (6004),  
Side rails (pair) (9002),  
Cushions for the side rails (pair) (6002P),  
Removable head/knee support (6007),  
Single footswitch for NV2560 (6200)  
Comes with a current Technical Passport issued,

### Warranty:

Lifetime warranty on frame,  
5 years on hydraulic and electric cylinders,  
5-year unconditional warranty on all parts, excluding upholstery damage,  
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE )



TAPICERKI STANDARD + (ZA DOPŁATK.)

