



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CLINNOVA Clinical 2 diagnostic and treatment table (NV2570-PRM-SEERSMEDICAL) (New)

New,
English production,
Very strong and rigid table available with a wide range of accessories,
Suitable for very heavy patients during examinations and minor medical procedures,
Height adjustment of the table and backrest by means of an electric actuator,
The table is equipped with a premium base and paper holder as standard,
Table advantages:

Retractable wheel design and adjustable leg for stability,
Frame construction designed for easy use of mobile lift,
Antimicrobial powder-coated base frame,
Easy Clean 70 cm wide upholstery as standard,
Choice of 15 standard upholstery colours
The table rises from minimum to maximum height in only 15 seconds,
Handheld remote control fitted as standard with this model,

Technical data:

Safe working load and lifting capacity of 265 kg,
Adjustable backrest from -20° to 80°,
Standard width of the bed is 70cm made of foam,
Table height adjustment range from 46 cm to 95 cm,

Additional accessories (at extra charge):

Face hole with plug in the bolster (6004),
Side rails (pair) (9002),
Cushions for the side rails (pair) (6002P),
Removable head/knee support (6007),
Double foot switch for NV2570 (6207),
Comes with a current Technical Passport issued,

Warranty:

Lifetime warranty on frame,
5 years on hydraulic and electric cylinders,
5-year unconditional warranty on all parts, excluding upholstery damage,
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATKĄ)

