



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Diagnostic and examination table CLINNOVA Mobile 2 (NV2450-CLS-SEERSMEDICAL) (New)

New,
English production,
Very strong and rigid table available with a wide range of accessories,
Suitable for very heavy patients during examinations and minor medical procedures,
Height adjustment of the table by a hydraulic cylinder, manual adjustment of the backrest by a gas spring,
The table is equipped with a classic base, side rails and paper holder as standard,
Table advantages:

Large castors with 125 mm individual brake (IBC),
Frame construction designed for easy use of the mobile elevator,
Base frame powder-coated with antimicrobial coating,
Easy Clean 70 cm wide upholstery as standard,
Choice of 15 standard upholstery colours

Technical data:


Safe working load and lifting capacity of 265 kg,
Adjustable backrest from -20° to 80°,
Standard width of the bed is 70 cm made of foam,
Table height adjustment range from 46 cm to 95 cm,

Additional accessories (at extra charge):

Face hole with plug in the bolster (6004),
Cushions for the side rails (pair) (6002P),
Removable head/knee support (6007),
Telescopic IV stand (6054),
Comes with a current Technical Passport issued,

Warranty:

Lifetime warranty on frame,
5 years on hydraulic and electric cylinders,
5-year unconditional warranty on all parts, excluding upholstery damage,
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE)

