



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **CLINNOVA diagnostic and examination table Mobile 2 (NV2660-CLS-SEERSMEDICAL) (New)**

New,  
English production,  
Very strong and rigid table available with a wide range of accessories,  
Suitable for very heavy patients during examinations and minor medical procedures,  
Height adjustment with electric actuator, manual backrest adjustment with gas spring,  
The table is equipped with a classic base, side rails and paper holder as standard,  
Table advantages:

Large castors with 125 mm central brake (LWMD),  
Frame construction designed for easy use of the mobile elevator,  
Base frame powder-coated with antimicrobial coating,  
Easy Clean 70 cm wide upholstery as standard,  
Choice of 15 standard upholstery colours,

Table rises from minimum to maximum height in just 15 seconds,

Technical data:

Safe working load and lifting capacity of 265 kg,  
Adjustable backrest from -20° to 80°,  
Standard width of the bed is 70 cm made of foam,  
Table height adjustment range from 46 cm to 95 cm,

Additional accessories (at extra charge):

Face hole with plug in the bolster (6004),  
Cushions for the side rails (pair) (6002P),  
Removable head/knee support (6007),  
Telescopic IV stand (6054),

Single footswitch for NV2460 and NV2660 (6200),

Comes with a current Technical Passport issued,

Warranty:

Lifetime warranty on frame,  
5 years on hydraulic and electric cylinders,  
5-year unconditional warranty on all parts, excluding upholstery damage,  
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD +[ZA DOPŁATK]

