



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### CLINNOVA Clinical 3 examination and treatment table (NV3560-PRM-SEERSMEDICAL) (New)

New,  
English production,  
Very strong and rigid table available with a wide range of accessories,  
Suitable for very heavy patients during examinations and minor medical procedures,  
Height adjustment by electric actuator, manual adjustment of backrest and footrest by gas spring,  
The table is equipped with a premium base,

#### Table advantages:

Retractable wheel design and adjustable leg for stability (RWD),  
Frame construction designed for easy use of the mobile lift,  
Base frame powder-coated with antimicrobial coating,  
Easy Clean 70 cm wide upholstery as standard,  
Choice of 15 standard upholstery colours,  
Table rises from minimum to maximum height in just 15 seconds,

#### Technical data:

Safe working load and lifting capacity of 265 kg,  
Adjustable backrest from -20° to 80° with anti-trap design,  
Adjustable footrest section 0° to 85° with anti-trap design,  
Standard width of the bed is 70 cm made of foam,  
Adjustable height of the table from 46cm to 95cm,

#### Additional accessories (extra charge):

Face hole with plug (6004),  
Side rails (pair) (9003),  
Cushions for the side rails (pair) (6003P),  
Upper limb stabiliser (6006),  
Removable head/knee support (6007),  
Single footswitch for NV3560 (6200),  
Comes with a current Technical Passport issued,

#### Warranty:

Lifetime warranty on frame,  
5 years on hydraulic and electric cylinders,  
5-year unconditional warranty on all parts, excluding upholstery damage,  
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD (+ZA DOPŁATEM)

