



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CLINNOVA Clinical 3 examination and treatment table (NV3570-CLS-SEERSMEDICAL) (New)

New,
English production,
Very strong and rigid table available with a wide range of accessories,
Suitable for very heavy patients during examinations and minor medical procedures,
Height, backrest and footrest adjustment by electric actuators,
The table is equipped with a classic base,

Table advantages:

Retractable wheel design and adjustable leg for stability (RWD),
Frame construction designed for easy use of mobile lift,
Base frame powder-coated with antimicrobial coating,
Easy Clean 70 cm wide upholstery as standard,
Choice of 15 standard upholstery colours,

Technical data:

Safe working load and lifting capacity of 265 kg,
Adjustable backrest from -20 ° to 80 °,
Adjustable footrest segment from 0° to 85°,
Standard width of the bed is 70cm made of foam,

Adjustable height of the table from 46 cm to 95 cm,

Additional accessories (at extra charge):

Face hole with plug (6004),
Side rails (pair) (9003),
Cushions for the side rails (pair) (6003P),
Upper limb stabiliser (6006),
Removable head/knee support (6007),
Triple footswitch for NV3570 (6209),
Comes with a current Technical Passport issued,

Warranty:

Lifetime warranty on frame,
5 years on hydraulic and electric cylinders,
5-year unconditional warranty on all parts, excluding upholstery damage,
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD (ZA DOPŁATKĄ)

