



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Diagnostic and examination table CLINNOVA Mobile 3 (NV3460-PRM-SEERSMEDICAL) (New)

New,
English production,
Very strong and rigid table available with a wide range of accessories,
Suitable for very heavy patients during examinations and minor medical procedures,
Height adjustment by electric actuator, manual adjustment of backrest and footrest by gas spring,
The table is equipped with a premium base, side rails and paper holder as standard,
Table advantages:

Table advantages:

- Large castors with 125 mm individual brake (IBC),
- Frame construction designed for easy use of the mobile elevator,
- Base frame powder-coated with antimicrobial coating,
- Easy Clean 70 cm wide upholstery as standard,
- Choice of 15 standard upholstery colours,
- Table rises from minimum to maximum height in just 15 seconds,

Technical data:

- Safe working load and lifting capacity of 265 kg,
- Adjustable backrest from -20° to 80° with anti-trap design,
- Adjustable footrest section 0° to 85° with anti-trap design,
- Standard width of the bed is 70 cm made of foam,
- Adjustable height of the table from 46cm to 95cm,

Additional accessories (extra charge):

- Face hole with plug (6004),
- Side rail pads (pair) (6003P),
- Upper limb stabiliser (6006),
- Removable head/knee support (6007),
- Telescopic IV-Pole (6054),

It has a current Technical Passport issued,

Warranty:

- Lifetime warranty on frame,
- 5 years on hydraulic and electric cylinders,
- 5-year unconditional warranty on all parts, excluding upholstery damage,
- Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁAT)

