



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynaecological chair CLINNOVA Gynae Pro (NV8583-CLS-SEERSMEDICAL) (New)

New,
English manufacture,
Adjustable seat height, backrest, footrests, tilt and CPR by electric actuator,
Retractable wheel design and adjustable leg for stability (RWD),
The chair comes with a classic base, 6041 side armrests, 6045 leg supports and NV9005 paper roll holder as standard,
Frame structure designed for easy use of mobile lift,
The table is equipped with a digital handheld remote control as standard,
Seat Benefits:
Electric height adjustment from 53 to 99 cm in just 15 seconds,
The chair has Trendelenburg and Anti Trendelenburg position (+/- 15°),
Remote control with digital display provides live information on height (cm), tilt angle (°) and Trendelenburg angle (°). Foot switch operation is available as an option (at extra cost),
Lifetime warranty on frame and 5 years warranty on all parts, excluding upholstery,
Choice of 15 standard antibacterial upholstery colours,
Antimicrobial powder-coated frame,
Free space under frame for portable lift,
Crescent shaped seat design that can be lifted for cleaning,
Fully adjustable patient leg rests,
Technical data:
70 cm bed width made of foam,
Safe working load and lifting capacity 285kg ,
Backrest adjustable from 0 ° to 80 °,
Adjustment range of table height from 53 cm to 99 cm,
Additional accessories (at extra charge):
Folding leg extension (NV9042)
Removable leg extension (NV9043)
Metal stirrups for gynaecological treatments (6044)
8-button footswitch (6210)
Comes with a current Technical Passport issued,
Warranty:
Lifetime warranty on frame,
5 years on hydraulic and electric cylinders,
5-year unconditional warranty on all parts, excluding upholstery damage,
Finance options (Poland only): Installments, Leasing, Loan

