



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CLINNOVA Podiatry Pro diagnostic and treatment chair (NV0495-PRM-SEERSMEDICAL) (New)

New,

English manufacture,

This chair has both an adjustable backrest with patient armrests and double adjustable footrests,

It is suitable for a wide range of seated examination procedures such as podiatry, lower limb ulcer clinics, orthopaedic casting (plaster surgery) and is also ideal for doctors' surgeries and health centres.

Adjustment of chair height, backrest, footrests, tilt and CPR by electric actuator,

The chair is equipped with large, individually braked IBC castors and a premium base,

The table is equipped with two adjustable armrests,

Frame structure designed for easy use of the mobile lift,

The foot section is adjustable and can be extended which is useful for taller patients,

The table comes with a handheld remote control as standard,

Benefits of the Chair:

Electric height adjustment from 53 to 99 cm in just 15 seconds,

The chair has Trendelenburg and Anti Trendelenburg positions (+/- 15°),

Folding armrests and NV9005 paper roll holder as standard,

Remote control with digital display provides live information on height (cm), tilt angle (°) and Trendelenburg angle (°). Foot switch operation is available as an option (at extra cost),

Lifetime warranty on frame and 5 years warranty on all parts, excluding upholstery,

Choice of 15 standard antibacterial upholstery colours,

Antimicrobial powder-coated frame,

Free space under the frame for a portable lift,

Technical data:

70 cm bed width made of foam,

Safe working load and lifting capacity 285kg ,

Backrest adjustable from 0° to 80° ,

Adjustment of footrest from 0 ° to -90 ° ,

Adjustment range of table height from 53 cm to 99 cm,

Additional accessories (at extra charge):

Fitted operator's chair (6100),

Fitted operator's premium chair (6150),

Battery back-up facility (6057)

Face hole with plug in the head rest (6004),

12 - button footswitch for CLINNOVA Podiatry Pro (NV9212),

Comes with a current Technical Passport issued,


Warranty:

Lifetime warranty on frame,

5 years on hydraulic and electric actuators,

5-year unconditional warranty on all parts, excluding upholstery damage,

Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD +[ZA DOPŁATK.]

