



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Tested disposable breathalyser TEST&DRIVE Sober Morning. 1 pc (New)**

New,  
Made in Poland,  
Packaging: cardboard box,  
TEST&DRIVE is a disposable alcotest for measuring the alcohol content in exhaled air,  
Attested,  
Effectiveness: 98 %,  
Reading time: 2 min,  
Easy to use,  
Allows accurate and fast measurement of the alcohol content in exhaled air,  
Professional, for use by every driver - as a supplement to the first aid kit in the car,  
Can also be effectively used by public institutions such as the police, municipal police, fire brigade, army, health services, etc.  
It is also perfect for all employers who want to make sure that their employees are safe, and that they do not endanger anyone or anything,  
Test&Drive consists of:  
Glass tube,  
Active substance in the form of adsorption granules,  
Filters,  
Aluminium foil,  
2 protective caps,  
User manual in languages: PL, GB, DE,  
Dimensions: 98 x 11 mm

