



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CLINNOVA Clinical 2 rehabilitation table (NV2550-CLS-SEERSMEDICAL) (New)

New,
English production,
Very strong and rigid table available with a wide range of accessories,
Suitable for very heavy patients during examinations and minor medical procedures,
Height adjustment of the table by a hydraulic cylinder, manual adjustment of the backrest by a gas spring,
The table is equipped with a classic base and paper holder as standard,
Table advantages:
Retractable wheel design and adjustable leg for stability,
Frame construction designed for easy use of mobile lift,
Antimicrobial powder-coated base frame,
Easy Clean 70 cm wide upholstery as standard,
Choice of 15 standard upholstery colours
Technical data:
Safe working load and lifting capacity of 265 kg,
Adjustable backrest from -20° to 80° with anti-trap design,
Standard width of the bed is 70 cm made of foam,
Table height adjustment range from 49 cm to 99 cm,
Additional accessories (at extra charge):
Face hole with plug in the bolster (6004),
Side rails (pair) (9002),
Cushions for the side rails (pair) (6002P),
Removable head/knee support (6007),
Comes with a current Technical Passport issued,
Warranty:
Lifetime warranty on frame,
5 years on hydraulic and electric cylinders,
5-year unconditional warranty on all parts, excluding upholstery damage,
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATĄ)

