



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)


Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## CLINNOVA Clinical 2 rehabilitation table (NV2550-PRM-SEERSMEDICAL) (New)

New,  
English production,  
Very strong and rigid table available with a wide range of accessories,  
Suitable for very heavy patients during examinations and minor medical procedures,  
Height adjustment of the table by a hydraulic cylinder, manual adjustment of the backrest by a gas spring,  
The table is equipped with a premium base and paper holder as standard,  
Table advantages:  
Retractable wheel design and adjustable leg for stability,  
Frame construction designed for easy use of mobile lift,  
Antimicrobial powder-coated base frame,  
Easy Clean 70 cm wide upholstery as standard,  
Choice of 15 standard upholstery colours  
Technical data:  
Safe working load and lifting capacity of 265 kg,  
Adjustable backrest from -20° to 80° with anti-trap design,  
Standard width of the bed is 70 cm made of foam,  
Table height adjustment range from 49 cm to 99 cm,  
Additional accessories (at extra charge):  
Face hole with plug in the bolster (6004),  
Side rails (pair) (9002),  
Cushions for the side rails (pair) (6002P),  
Removable head/knee support (6007),  
Comes with a current Technical Passport issued,  
Warranty:  
Lifetime warranty on frame,  
5 years on hydraulic and electric cylinders,  
5-year unconditional warranty on all parts, excluding upholstery damage,  
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE )

