



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Rehabilitation table CLINNOVA Mobile 3 (NV3470-CLS-SEERSMEDICAL) (New)

New,  
English production,  
Very strong and rigid table available with a wide range of accessories,  
Suitable for very heavy patients during examinations and minor medical procedures,  
Height, backrest and footrest adjustment by electric actuator,  
The table is equipped with a classic base, side rails and paper holder as standard,  
Table advantages:  
Large castors with 125 mm individual brake (IBC),  
Frame construction designed for easy use of the mobile elevator,  
Base frame powder-coated with antimicrobial coating,  
Easy Clean 70 cm wide upholstery as standard,  
Choice of 15 standard upholstery colours,  
Table rises from minimum to maximum height in just 15 seconds,  
Technical data:  
Safe working load and lifting capacity of 265 kg,  
Adjustable backrest from -20 ° to 80 °,  
Adjustable footrest segment from 0° to 85°,  
Standard width of the bed is 70cm made of foam,  
Adjustable height of the table from 46 cm to 95 cm,  
Additional accessories (at extra charge):  
Face hole with plug (6004),  
Cushions for the side rails (pair) (6003P),  
Upper limb stabiliser (6006),  
Removable head/knee support (6007),  
Telescopic IV-Pole (6054),  
It has a current Technical Passport issued,  
Warranty:  
Lifetime warranty on frame,  
5 years on hydraulic and electric cylinders,  
5-year unconditional warranty on all parts, excluding upholstery damage,  
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATĄ)

