



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Rehabilitation table CLINNOVA Mobile 3 (NV3670-CLS-SEERSMEDICAL) (New)

New,  
English production,  
Very strong and rigid table available with a wide range of accessories,  
Suitable for very heavy patients during examinations and minor medical procedures,  
Height, backrest and footrest adjustment by electric actuator,  
The table is equipped with a classic base, side rails and paper holder as standard,  
Table advantages:

Large castors with 125 mm central brake (LWMD),  
Frame construction designed for easy use of the mobile elevator,  
Base frame powder-coated with antimicrobial coating,  
Easy Clean 70 cm wide upholstery as standard,  
Choice of 15 standard upholstery colours,  
Table rises from minimum to maximum height in just 15 seconds,

### Technical data:

Safe working load and lifting capacity of 265 kg,  
Adjustable backrest from -20 ° to 80 °,  
Adjustable footrest segment from 0° to 85°,  
Standard width of the bed is 70cm made of foam,  
Adjustable height of the table from 46 cm to 95 cm,

### Additional accessories (at extra charge):

Face hole with plug (6004),  
Cushions for the side rails (pair) (6003P),  
Upper limb stabiliser (6006),  
Removable head/knee support (6007),  
Telescopic IV-Pole (6054),

It has a current Technical Passport issued,

### Warranty:

Lifetime warranty on frame,  
5 years on hydraulic and electric cylinders,  
5-year unconditional warranty on all parts, excluding upholstery damage,  
Finance options (Poland only): Installments, Leasing, Loan



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATKĄ)

