



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



---

### **Docking station for 6 pumps FRESENIUS VIAL ORCHESTRA BASE A (Reconditioned)**

---

Used,  
Technical condition: very good,  
Visual condition: very good,  
French production,  
Power supply: 230 V,  
Emergency power supply: new built-in 6 V 1.3 Ah battery  
Frequency: 50/60 Hz,  
Power: 125 VA,  
The Base A docking station is a data hub and power supply for up to 6 DPS modules and/or DPS Visio modules and/or MVP modules,  
Allows faster therapeutic control,  
Communication interface compatible with the main PDMS software,  
Built-in RS232 interface,  
One single power cable for all connected modules,  
Minimum battery life of 2 hours,  
Equipped with tripod with comfortable carrying handle,  
No pumps included - can be retrofitted to FRESENIUS ORCHESTRA infusion pumps (syringe or drip) at additional cost  
Operating instructions in German (PDF),  
Docking station dimensions:  
Height: 8 cm,  
Width: 34.5 cm,  
Depth: 25.5 mm,  
Height including tripod: 80 cm,  
Weight of station with tripod: 4 kg,  
Weight of docking station with pumps: 17.5 kg,  
It has an up-to-date Technical Passport issued,  
Warranty:  
6 months for domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan

