



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment table Wesseling BASIX 2 (BAS2TEMW) (New)

New,
Dutch manufacture,
This two-piece treatment table is suitable for various purposes.
Thanks to the scissor frame, the table has a wide range of height adjustment and lifts straight up.
The headrest can be adjusted both positively and negatively.
As standard, the table is adjustable with an electric actuator, has wheels with braking system and Trendelenburg position,
Technical parameters:
Length: 200 cm,
Width: 64 cm,
Upholstery thickness: 50 mm
Height adjustment: from 45 to 115 cm,
Headrest adjustment: -35° to 85°,
Adjustment of reclining section: 0° to -25°,
Load capacity : 225 kg,
Safe working load: 250 kg,
Duration of height adjustment: 40 sec.
Adjustable sections: 2,
Brake system: central,
Wheels: Ø 75 mm,
As additional options available (at extra charge):
Paper roll holder,
Nose hole,
Large neck cushion,
Foot height adjustment (control from the frame),
Special table top width 70cm,
Special top width 80cm,
Thicker upholstery 70mm,
Terry cloth cover for tables 65-70 cm,
Terry cloth cover for tables 65-70cm with cut-out for face.
It has an up-to-date Technical Passport,
Warranty: 24 months

