



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment table Wesseling BASIX 3 FLEX (BAS3FHMW) (New)

New,

Dutch manufacture,

This three-piece treatment table is suitable for a variety of purposes.

Equipped with a nose hole and a rod switch system, the table is the ideal table for physiotherapy,

Due to the division of the table surface, the patient always sits or lies in a comfortable position during treatment.

In addition, the highest position of 115cm ensures that everyone can create an ergonomic working surface

The table is equipped with two gas springs under the central section, which means that you can achieve the so-called "kidney upright position",

As standard, the table is adjustable with a hydraulic cylinder, has wheels with braking system and drainage and Trendelenburg positions,

Technical parameters:

Length: 200 cm,

Width: 64 cm,

Upholstery thickness: 50 mm

Height adjustment: 48 to 98 cm,

Adjustment of the section of the bed: 0° to -25°,

Headrest adjustment: -55° to 85°,

Backrest adjustment: -25° to 85°,

Load capacity : 225 kg,

Safe working load: 250 kg,

Duration of height adjustment: 40 sec.

Adjustable sections: 3,

Brake system: central,

Wheels: Ø 75 mm,

As additional options available (at extra charge):

Paper roll holder,

Nose hole,

Large neck cushion,

Foot height adjustment (control from the frame),

Special table top width 70cm,

Special top width 80cm,

Thicker upholstery 70mm,

Terry cloth cover for tables 65-70 cm,

Terry cloth cover for tables 65-70cm with cut-out for face.

It has an up-to-date Technical Passport,

Warranty: 24 months

